

ODM Full Arch Zirconia Implants Bridge Zirconium Implants Holistic Dentistry

Our Product Introduction

for more products please visit us on adschinadentalallab.com


Basic Information

- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Dental box and poly bag
- Delivery Time: 3-4 Days
- Payment Terms: T/T, PayPal, Western Union
- Supply Ability: 200 Technicians



Product Specification

- Abutment Height: 3mm
- Implant Length: 10mm
- Abutment Type: Screw-Retained
- Thread Pitch: 1.5mm
- Platform Height: 2.5mm
- Surface Treatment: Sandblasted And Acid-Etched
- Connection Type: Internal Hex
- Material: Titanium
- Abutment Angle: 15 Degrees
- Implant Diameter: 3.5mm
- Thread Angle: 60 Degrees
- Platform Diameter: 3.75mm
- Thread Type: Self-Tapping
- Thread Depth: 1.2mm
- Highlight: ODM full arch zirconia implants .



More Images



Product Description

Advantages

Our

- Secure and reliable restoration
- Optimal function and esthetics
- Customized fit for optimal function and esthetics
- Versatility for various implant restoration scenarios
- Collaboration with dental professionals for seamless integration

Key Features

****Secure and Reliable:**** Screw-Retained Implants utilize a screw attachment to secure the implant restoration in place. The screw provides a secure and reliable connection between the implant abutment and the prosthetic restoration, ensuring stability and longevity.

****Optimal Function and Esthetics:**** Our Screw-Retained Implant restorations are designed to optimize both function and esthetics. They are customized to blend seamlessly with the patient's natural teeth, providing a harmonious and natural-looking smile while ensuring proper occlusion and chewing efficiency.

****Customized Fit:**** Each Screw-Retained Implant restoration is precisely fabricated to fit the unique anatomy of the patient's mouth. Our experienced technicians work closely with dental professionals to ensure a customized fit that promotes optimal function and esthetics.

****Versatility:**** Screw-Retained Implants can be used for single-tooth replacements as well as for multiple teeth, including implant-supported bridges. They offer versatility in restoring both anterior and posterior regions of the mouth.

****Collaboration with Dental Professionals:**** We maintain open communication and collaborate closely with dentists and oral surgeons to ensure the seamless integration of the Screw-Retained Implant restoration within the overall treatment plan. This collaborative approach ensures optimal patient outcomes.

Indications

- Single-tooth replacement
- Multiple teeth replacement (including implant-supported bridges)

Considerations

- Proper oral hygiene and regular dental visits are essential for the long-term success and maintenance of Screw-Retained Implant restorations.
- Screw-Retained Implants may require occasional access for maintenance and adjustments, which can be done easily with the removal and reattachment of the restoration.